

Pick Up **New** Products

Sanplatec Antistatic Product Series

Product name:
Antistatic Bottle, Antistatic PE dispo cup,
Antistatic round washed bottle



■ Features

We added new function Antistatic to normal PE bottle.

Special polymer antistatic agent allows semi-permanent antistatic effect

The antistatic agent is contained to main body, inner seal and cap.

(Capability body→surface resistivity (Ω) : $10^8 \sim 10^{11}$)

(JIS K 6911 thermosetting plastic general test (5.13 resistivity) Only for reference)

■ Solve concerns regarding the bottle static electricity

- When handling powdered material, it sticks to inner wall.
- Static electricity affect the content
- Dust sticks easily

■ Caution

- The effect of this item is 「antistatic」 , not 「Static elimination」 .
- Do not use with organic solvent.
Please be careful to use acid and alkali which may cause performance degradation.