

Date: 10/04/23

MODEL 2501 XRST PORTASPEC® WDXRF SPECTROGRAPH

RJM Systems (Pty) Ltd,
PO Box 31347
Kyalami, 1684 South Africa
251 Papenfus Drive, Beaulieu, Kyalami.
Midrand 1684



WHAT IS INCLUDED



- 2522X Water Coolant Circulator, complete with Pump, Hoses, Fittings, and Coolant Mixture.
- PC Interface with internal mini-computer hard drive and PortaSpec “PCS” Operating Software.
- Computer Monitor and Mouse/Keyboard.
- Catalog No. 2526X Premium Shipping Container, molded plastic with custom foam-lined insulation larger in size to provide extra protection for shipping. Accommodates PortaSpec Instrument **only**.

PortaSpec Series X WDXRF Spectrograph

Specifications:

- 115 Volts, 60 Hertz, Single Phase, 30 kVCP, 5 mA.
- Complete with integrated Digital Scaler and Timer, Tungsten Target X-ray Tube
- Standard Proportional Detector, LiF 200 Crystal, mounted on a rotating Goniometer with the ability to pre-set up to five (5) elements of interest with two (2) set for Cr (Chrome) and Ti (Titanium)
- Rotating Sample Tray that can hold up to eight (8) samples, including Blank/Standards, for sequential analysis, up to a maximum sample size of 2-9/16" dia. x 1" thick.

OPTIONAL ACCESSORY (STRONGLY RECOMMENDED):

MCR SOLA POWER LINE CONDITIONER

- To provide excellent noise filtering and surge suppression to protect PortaSpec from damage, degradation or mis-operation.
- Improves PortaSpec stability.
- Combined with excellent voltage regulation, the MCR can increase the actual “Mean Time Before Failure” (MTBF) of the PortaSpec.
- The MCR is a perfect choice where dirty power, caused by impulses, swells, sags, brownouts, and waveform distortion that can lead to costly downtime because of damage to the PortaSpec.
- 500 VA, 60 Hz, Input +10%/-20% Output: +/- 5

WARRANTY

- Two (2) Year Limited Warranty on new PortaSpec Equipment
- **Five (5) Year Pro-Rated Warranty on new X-ray Ray Tubes**
- Computer **Accessories** carry 1-year manufacturer's warranty only
- Some exclusions apply